



DNA from Saliva

EASY • PAINLESS • PROVEN

Collect high quality DNA samples for molecular diagnostics with an FDA cleared device[†].

Genetic testing starts with the collection of DNA samples. Use Oragene®Dx for proven collection, stabilization and transportation of DNA from saliva.

- **High quality, high quantity DNA[†]**
- **Fast and reliable self collection**
- **Long term room temperature sample stability**
- **Transport and store at room temperature**
- **Integrated 1D barcode + optional 2D bottom barcode for laboratory workflow efficiency**

Molecular diagnostic test manufacturers could reduce the time and regulatory burden of their own FDA submission by leveraging the 510(k) clearance of Oragene•Dx.



To request free sample kits, go to www.dnagenotek.com

For In Vitro Diagnostic Use
CE
 For collection of human DNA



DNA Genotek Inc.
 3000 - 500 Palladium Drive
 Ottawa, ON, Canada K2V 1C2
 Subsidiary of OraSure Technologies, Inc.

Toll-free (North America): 1.866.813.6354
 Tel.: +1.613.723.5757 • Fax: +1.613.723.5057
www.dnagenotek.com
info@dnagenotek.com

DNAgenotek


Your kit, your way.

The Oragene-Dx 610 combines ease-of-use and high quality samples with the ability to tailor the device to your needs. Tube length, anti-rotation features for automation and enhanced barcoding options give you ultimate flexibility for storage, sample processing and traceability. GenoFIND™ services are available to supplement your sample management needs.

Sample size	Recommended configuration	Laboratory advantages	Collection, transport and storage advantages	Supplemental services
Small to mid-size collection 	Standard 610 kit with 1D tube barcode	<ul style="list-style-type: none"> Compatible with manual processing Automated de-capping (low-medium throughput) Sample tracking with 1D tube barcode 	<ul style="list-style-type: none"> Non-invasive At home or centralized collection Quick and easy self collection Instructions available in multiple languages Ship via standard postal system No cold chain required Compact design minimizes storage space 	<ul style="list-style-type: none"> Custom sample packaging for brand awareness Complete fulfillment capability Regulatory guidance Mailing accessories Wet lab services in a GLP/CLIA accredited lab
Large-scale collection 	Standard 610 kit with 1D tube barcode • Optional matching 2D bottom barcode*	<ul style="list-style-type: none"> Automated de-capping (high-throughput) Sample tracking with option for bed scanning 		

* No false bottom on configurations with 2D bottom barcodes

Kit contents and packaging specifications

For full specifications on weight and dimensions, refer to PD-BR-00233 on our website.

Dimensions: 14.0 x 7.5 x 2.8 cm

Weight: 37 g

Shelf-life: 30 months at room temperature

Sample stability: 12 months



Collector



Instructions



Tube cap



Packaging

†Performance characteristics (representative data)**

	Yield (µg)	Concentration (ng/µL)	A _{260/280} ratio
Mean ± SD	36.0 ± 41.1	76.3 ± 78.0	1.7 ± 0.1
Median	27.0	52.1	1.7
90% of samples	≥ 7.5	≥ 16.0	≥ 1.6

** Data obtained by testing 45 samples from 45 individuals.

Oragene-Dx has been cleared for in vitro diagnostic use in the U.S.A.






† FDA cleared for in vitro diagnostic use; for prescription and direct-to-consumer use.

*Oragene is a registered trademark, and GenoFIND is a trademark of DNA Genotek Inc.

All other brands and names contained herein are the property of their respective owners.

Patent (www.dnagenotek.com/legalnotices)

Easy collection



Full collection instructions can be found on our website.

GENOFIND™
Services
Superior support, analysis and reporting