

**Instructions for sample collection:**



**Collection precautions:**

**Ensure the sponge tip does NOT come into contact with any surface prior to collection.**

**Donor should NOT eat, drink, smoke or chew gum for 30 minutes before collecting saliva sample.**

**Intended use:** For the collection and stabilization of RNA from human saliva samples.

**Contents:** Contains 1 collection kit.

**Warning and precautions:** Choking hazard. Caution should be used when inserting sponge into the mouth.

Wash with water if liquid comes in contact with eyes or skin. Do NOT ingest. See MSDS at [www.dnagenotek.com](http://www.dnagenotek.com)

**Storage:** 15-25°C

**Summary and explanation of the kit:**

ORACollect®-RNA is a self-collection kit that provides the materials and instructions for collecting human saliva samples.

**Label legend:**

- Catalog number
- CE Marking
- In vitro diagnostic medical device
- Manufacturer
- Authorized Representative
- 15°C-25°C Storage instructions
- Collect sample by (Use by)
- Lot number
- For Prescription Use Only

	<p>Open package and remove collector without touching sponge tip. Place sponge as far back in the mouth as comfortable and rub along the lower gums (see close up image) in a back and forth motion. Gently rub the gums 10 times. If possible, avoid rubbing the teeth.</p>
	<p>Gently repeat rubbing motion on the opposite side of the mouth along the lower gums for an additional 10 times.</p>
	<p>Hold the tube upright to prevent the liquid inside the tube from spilling. Unscrew the cap from the collection tube without touching the sponge.</p>
	<p>Turn the cap upside down, insert the sponge into the tube and close cap tightly.</p>
	<p>Invert the capped tube and shake vigorously 15 times.</p>



For IVD, Rx and For Use Under Emergency Use Authorization.

This sample collection device has not been FDA cleared or approved.

This sample collection device has been authorized by FDA under an EUA.

This sample collection device has been authorized only for the collection and maintenance of saliva specimens as an aid in detection of nucleic acid from SARS-CoV-2, not for any other viruses or pathogens.

This sample collection device is only authorized for the duration of the declaration that circumstances exist justifying the authorization of emergency use of in vitro diagnostics for detection and/or diagnosis of COVID-19 under Section 564(b)(1) of the Federal Food, Drug and Cosmetic Act, 21 U.S.C. § 360bbb-3(b)(1), unless the authorization is terminated or revoked sooner.



*Superior samples • Proven performance*

■ DNA Genotek Inc.  
3000 - 500 Palladium Drive  
Ottawa, ON, Canada K2V 1C2

Toll-free (North America): 1.866.813.6354  
Tel: +1.613.723.5757 • Fax: +1.613.723.5057  
info@dnagenotek.com • www.dnagenotek.com



Emergo Europe, Prinsessegracht 20, 2514 AP The Hague, The Netherlands

Australian Sponsor: Emergo Australia, Level 20, Tower II, Darling Park, 201 Sussex Street, Sydney, NSW 2000 Australia

\*ORACollect is a registered trademark and DNA Genotek is a trademark of DNA Genotek Inc. Some DNA Genotek products may not be available in all geographic regions, contact your sales representative for details. All DNA Genotek protocols, white papers and application notes, are available in the support section of our website at [www.dnagenotek.com](http://www.dnagenotek.com).

Patent ([www.dnagenotek.com/legalnotices](http://www.dnagenotek.com/legalnotices))  
© 2020 DNA Genotek Inc., a subsidiary of OraSure Technologies, Inc., all rights reserved.  
PD-PR-00764 Issue 4/2020-10