

Isifinyezo nencazelo yekhithi:

OMNIGENE-GUT inikezela ngezinto nemiyalo yokuqoqa nokuqinisa i-microbial DNA kusuka kusampuli ye-fecal.

Izexwayiso nokuqaphela:

- NGOKUSEBENZISA KWANGAPHANDLE KUPHELA.
- UNGALOKOTHI ususe okuphezulu kweshubhu okuphuzi kusuka kushubhu.
- UNGALOKOTHI ucithe amanzi okuqinisa eshubhini.
- Hlanza ngamanzi uma okusamanzi kuhlangana namehlo noma isikhumba sakho. UNGALOKOTHI ugwinye.
- Uma uqoqa isampuli ye-fecal yokusamanzi, bona imiyalo yomsebenzisi enikezelwe ngokwehlukile.
- Izinto ezincane zingabangela ubungozi bokuvaleka.

Isitoreji: 15°C to 25°C

Hambisa ngokuya ngemithetho ebekiwe emboza isithuthi samasampuli ofuzo. Bona i-MSDS ku-www.dnagenotek.com

Ilejendi yelebuli:

- Qoqa isampuli ngo- (Sebenzisa ngo-)
- Inombolo yekhathalogi
- Umkhiziqi
- 15°C to 25°C Imiyalo yesitoreji
- Isexeyiso, xhumana nemiyalo yokusebenzisa
- Inombolo ye-Lot

IMIYALO YOMSEBENZISI

Funda yonke imiyalo ngaphambi kokuqoqa

Inqubo:



- 1 IZIQAPHELISO EZIBALULEKILE:**
- Citha i-bladder yakho ngaphambi kokuqala ukuqoqa.
 - Qoqa esampuli ye-fecal engenawo umchamo noma amanzi wasethoyilethi.
 - I-Toilet paper noma ama-tissue angadingeka.




- 2** Ngenkathi ubambe okuphezulu kweshubhu okuphuzi, khumula KUPHELA isivalo esiphephuli kusuka kukhithi usibekele eceleni ukuze usebenzise kamuva.



OKUBALULEKILE:


UNGALOKOTHI ususe okuphezulu kweshubhu okuphuzi.
UNGALOKOTHI ucithe amanzi okuqinisa eshubhini.



- 3** Sebenzisa ispashula ukuqoqa inani lesampuli le-fecal.




Usayizi ovamile wesampuli ye-fecal.



- 4** Dlulisela isampuli ye-fecal phezu kushubhu okuphuzi. Phinda isampuli ize igcwalise okuphezulu kweshubhu okuphuzi.




OKUBALULEKILE: UNGALOKOTHI uphushela isampuli ngaphakathi kweshubhu.




- 5** Beka ngokuvundlile phzulu kushubhu ukuze ulinganise isampuli futhi ususe okunye okuningi.


Sula ingaphandle leshubhu nokuphezulu nge-toilet paper noma i-tissue njengokudingekayo.



- 6** Thatha isivalo esiphephuli esinomphetho oqinile ubheke phansi bese uvale okuphezulu kweshubhu okuphuzi kuze kuvaleke kuqine.



- 7** Xukuza ishubhu elivaliwe ngokuqina nangokushesha ngokuya phambili nasemuva ubuncane obungamasekhondi angu-30.



- 8** Isampuli ye-fecal izoxhutshwa nokusamanzi kokuqinisa eshubhini; akuwona wonke amacezu azoncibilika.

OKUBALULEKILE: Qhubeka nokuxukuza uma amacezu **amakhulu** asala njengoba kuboniswe esithombeni A.



- 9** Beka ispashula ephakhejini loqobo noma sivale nge-toilet paper bese ulahle ezibini.

OKUBALULEKILE: Thumela isampuli ukuze icutshungulwe ulandele imiyalo yokuletha enikezwe ngokuhlukile umhlinzeki wekhithi.