

Isishwankathelo kunye nenkcazelo yekiti:

OMNIGENE-GUT Unika izixhobo kunye nemiyalelo yokuqokelela nokuzinzisa i-microbial DNA esuka kiwsampula ye-fecal.

Izilumkiso kunye nokuqaphela:

- UKUSETYENZISWA NGAPHAKATHI KUNYE.
- UNGAKE ususe the ityhubhu elilubheli ngasephezulu kwetyhubhu.
- UNGAKE uchithe amanzi zinzisayo kwityhubhu.
- Hlamba ngamanzi ukuba amanzi aqhakamshelana nelihlo okanye isikhumba. UNGAKE uyitye.
- Ukuba uqokelela isampula ye-fecal, bona eyahlukileyo imiyalelo yomsebenzisi.
- Izinto ezincinane zingabangela ubungozi bokuxhelwa.

Ukugcinwa: 15°C to 25°C

Thumela ngokuhambelana nemigaqo evumelekileyo equka ukuthuthwa kwamachiza ezempilo yendayo. Bona i-MSDS e www.dnagenotek.com


Umbalo weqala:

- ☒ Qokelela isampula nge (Isetyenziswa ngu)
- REF Inombolo ye-Catalog
- 🏠 Umenzi
- 15°C / 25°C Imiyalelo yocino
- ⚠️ Qhaphela, qhakamshelana nemiyalelo yokusetyenziswa
- LOT Inombolo yeLot

USEIMIYALELO YOMSEBENZISI

Funda yonke imiyalelo ngaphambi kokuqokelela

Inkqubo:



- 1 AMALINGISELELO ABALULEKILEYO:**
- Citha i-bladder yakho ngaphambi koqokelelo.
 - Qokelela isampula ye-fecal mahala yomchamo okanye amanzi angasese.
 - Iphepha langasese okanye iithishu zingadingeka.




- 2** Ngelixa ubambe ityhubhu elilubheli phezulu, khulula KUPHELA isiciko esifusi ukusuka kwikithi kwaye uyibeke bucala uze uyisebenzise kamva.



KUBALULEKILE:


UNGAKE ususe ityhubhu elubheli phezulu. UNGAKE uchithe amanzi okuzinzisa kwityhubhu.



- 3** Sebenzisa i-spatula ukuqokelela umyinge omncinane wesampula ye-fecal.



Umlinganiso woqobo wesampula ye-fecal.




- 4** Tshintshela isampula ye-fecal kwityhubhu elubheli ephezulu. Phunda ide isampula icwalise ityhubhu elubheli phezulu.



KUBALULEKILE: UNGAKE ududule isampula kwityhubhu.



- 5** Krwela ngokulele kumphezulu wetyhubhu ukulinganisa isampula kwaye ususe eninzi. Sula ityhubhu ngaphandle nangaphezulu ngephepha langasese okanye ithishu xakudingeka.




- 6** Chola isiciko esimfusi esinesiphelo esiqinileo sijonge phantsi kwaye isibeke kwityhubhu elilubheli phezulu side sivalwe ngokuqinileyo.



- 7** Shukumisa ityhubhu ezaliweyo ngokuqina nokukhawuleza njengoko unokukwazi ngokushukuma kwasemva naphambili kangangobucinane b30.

30 mizuzwana.



- 8** Isampula ye-fecal izakuxutyhwa namanzi okuzinzisa kwityhubhu; akuzizo zonke iizinro ezizakunyibilika.

KUBALULEKILE: Qhubeka ushukumise ukuba ekulu iizinto zihlala njengoko ziboniswa ku Mzobo A.



- 9** Beka i-spatula kwigxowa yayo okanye uyijije ngephepha langasese kwaye uyilahle emgqomeni.

KUBALULEKILE: Thumela isampula iyohlolwa kulandela imiyalelo yothunyelwe eziswe ngokwahlukileyo ngunikezeli we-kit.