



Non-invasive, assisted collection of high quantity and quality genomic DNA from saliva of young children

J.O. Niles, A. Jackson, S. Rabuka, and R.M. Iwasiov
DNA Genotek, Ottawa, Ontario, Canada

Abstract

Large population-based studies, involving thousands of subjects, are increasingly being used to investigate the genetic determinants of complex diseases. Saliva is a convenient source of genomic DNA because it can be collected in a painless and non-invasive manner. Non-invasive methods and techniques that permit self-collection are preferred because they increase compliance rates and reduce costs. For this reason, many large-scale studies now use DNA from saliva collected using the Oragene® kit† as the source of genomic DNA for downstream applications. Human genomic DNA from saliva can be used for genotyping, sequencing and micro-array analysis. In order to facilitate the non-invasive collection of genomic DNA from a population of all ages (including those individuals who cannot spit) we describe a new method of using sponges to transfer saliva from a donor's mouth into an Oragene kit. Unlike buccal swabs, this method transfers saliva which contains high quality and quantity DNA into the Oragene kit which preserves the DNA and prevents bacterial growth. We report that this method allows for the collection of a median yield of 17.3 µg of genomic DNA with median A_{260}/A_{280} ratio of 1.8 and a molecular weight > 23 kb. The collected samples are stable at ambient temperature for years and can be used in a multiple downstream applications.

Materials and methods





	1 Place one sponge in cheek pouch. Gently move the sponge along the gums and inner cheeks for 30 seconds to soak up as much saliva as possible.
	2 Once the sponge is saturated with saliva, insert sponge in V-notch of funnel. Wring saliva out of sponge using a twisting and pushing motion against the inner wall of the V-notch. Saliva will flow into tube.
	3 Repeat these steps (1 to 2) USING THE SAME SPONGE until the liquid saliva (not bubbles) reaches the fill line. Check sponge for damage each time before inserting into donor's mouth. Use second sponge if first sponge shows any signs of wear or tear. Tap tube bottom against hard surface to reduce number of bubbles.
	4 Hold the tube upright with one hand. Close the lid with the other hand (as shown) by firmly pushing the lid until you hear a loud click. The liquid in the lid will be released into the tube to mix with the saliva. Make sure the lid is closed tightly.
	5 Hold the tube upright. Unscrew the funnel from the tube.
	6 Use the small cap to close the tube tightly.
	7 Shake the capped tube for 5 seconds. Discard or recycle the funnel. Discard sponges.

Figure 1: Assisted collection of saliva using a sponge with the Oragene kit.

† Saliva samples were collected with Oragene®•DNA or Oragene®•DISCOVER.

Results




Figure 2: Distribution of donor ages

- Children aged 9 months to 10 years old were recruited for this study.
- Median age of children was 3.5 years.
- Samples from a total of 77 children were analyzed.




Figure 3: Time required to collect saliva to the fill-to-line

- 65% of all donors irrespective of age were able to complete the collection of saliva in under 6 minutes.




Figure 4: Total amount of DNA recovered.




Figure 5: Concentration of DNA purified from collected saliva.

- Collected samples were purified according to the prepIT™•L2P (DNA Genotek) purification protocol¹.
- Purified DNA was quantified by fluorescence using SYBR® Green I dye².
- The total amount of DNA collected from each child is reported in Figure 4. The median amount of DNA recovered was 17.3 µg.
- The resulting purified DNA had a median concentration of 60.2 ng/µL.




Figure 6: Corrected A₂₆₀/A₂₈₀ ratio of purified DNA.

- The corrected A₂₆₀/A₂₈₀ ratio was calculated by subtracting the A₃₂₀ value (which represents scattered light due to insoluble material) from both the A₂₆₀ and A₂₈₀ values.
- The purified DNA had a median A₃₂₀ corrected A₂₆₀/A₂₈₀ ratio of 1.82.




Figure 7: Representative agarose gel showing high molecular DNA.

- The molecular weight of purified DNA was assessed by running a 0.8% agarose gel. The purified DNA consistently had a molecular weight > 23 kb as compared with the Lamda-Hind III ladder.

