

A real-time PCR assay was used to demonstrate stability of human RNA in Oragene•RNA/saliva samples stored at room temperature (Figure 2). The data shows that the Oragene•RNA solution stabilized RNA messages for 8 weeks as demonstrated by the low variability in Ct values. Real-time PCR signal for each of the 5 specific RNA messages did not deviate more than 2 Ct values (Figure 3). Results summarize the data for 5 gene products of interest: human 18S rRNA, β -actin, β -2-microglobulin, Interleukin-8, and Histatin-3.




Figure 2: Real-time PCR data for 5 Oragene•RNA/saliva samples stored at room temperature for 1 (solid lines) and 8 weeks (dashed lines) before being purified and probed using primers for human 18S rRNA. RNA from cultured HCT116 cells was used as a positive control.

Donor	Time (weeks)	Ct value for human 18S rRNA	-RT Ct value for human 18S rRNA	Ct value for human β -actin	-RT Ct value for human β -actin
A	1	12.6	26.7	23.6	29.8
	8	13.0	29.1	23.2	30.2
B	1	13.7	27.9	23.2	32.6
	8	16.5	28.6	27.0	29.6
C	1	13.9	26.5	22.5	30.6
	8	13.8	27.8	23.0	31.9
D	1	13.7	28.4	26.2	31.8
	8	13.3	28.5	25.3	30.1
E	1	14.9	28.6	25.6	30.8
	8	13.6	29.1	22.7	28.1

Table 2: Summary of Ct values from real-time PCR data for human 18S rRNA and β -actin mRNA.




Figure 3: Summary of Ct values for 5 genes. RNA from saliva from 8 donors collected and stored up to 8 weeks was purified at given time points. Real-time PCR was used to detect the 5 genes of interest. Ct values at each time point were compared to Ct values obtained from samples purified immediately after collection. The red curve represents the average Ct for all donors.

Conclusions

Oragene•RNA is a non-invasive RNA self-collection kit that can be used by untrained subjects to collect high quality RNA. The RNA collected and purified from saliva using Oragene•RNA is high quality, with an average corrected A_{260}/A_{280} of 1.9. Oragene•RNA can stabilize total RNA at room temperature for up to 8 weeks as demonstrated by agarose gel electrophoresis. Specifically, 5 human specific genes were demonstrated by real-time PCR to be stable at room temperature for 8 weeks.

References

- 1 Oragene•RNA (RE-100) collection instructions. DNA Genotek. PD-PR-020.
- 2 Oragene•RNA purification protocol using the Qiagen RNeasy Micro Kit for volumes up to 1,000 μ L. DNA Genotek. PD-PR-021.

*Oragene is a registered trademark of DNA Genotek Inc. All other brands and names contained herein are the property of their respective owners. All DNA Genotek protocols, white papers and application notes, are available in the support section of our website at www.dnagenotek.com.