REF ORE-100

**DNagenotek™**

■ DNA Genotek Inc.  
3000 - 500 Palladium Drive  
Ottawa, ON, Canada K2V 1C2  
Tel.: +1.613.723.5757  
info@dnagenotek.com  
www.dnagenotek.com

**EN** ENGLISH

**Read all instructions prior to collection.**

**Collection precautions:**

Ensure the sponge tip does NOT come into contact with any surface prior to collection.

Donor should NOT eat, drink, smoke or chew gum for 30 minutes before collecting saliva sample.

**Intended use:** For the collection and stabilization of RNA from human saliva samples.

**Contents:** Contains 1 collection kit.

**Warning and precautions:**

Choking hazard. Caution should be used when inserting sponge into the mouth.

Wash with water if liquid comes in contact with eyes or skin. Do NOT ingest. See MSDS at www.dnagenotek.com

**Storage:** 15-25°C

**Summary and explanation of the kit:**

ORACollect-RNA is a self-collection kit that provides the materials and instructions for collecting human saliva samples.

**Instructions for sample collection:**

- 1 Open package and remove collector without touching sponge tip. Place sponge as far back in the mouth as comfortable and rub along the lower gums (see close up image) in a back and forth motion. Gently rub the gums 10 times. If possible, avoid rubbing the teeth.
- 2 Gently repeat rubbing motion on the opposite side of the mouth along the lower gums for an additional 10 times.
- 3 Hold the tube upright to prevent the liquid inside the tube from spilling. Unscrew the cap from the collection tube without touching the sponge.
- 4 Turn the cap upside down, insert the sponge into the tube and close cap tightly.
- 5 Invert the capped tube and shake vigorously 15 times.

**Label legend:**

- Catalog number
- CE Marking
- In vitro diagnostic medical device
- Manufacturer
- Authorized Representative
- Storage instructions
- Collect sample by (Use by)
- Lot number

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OraSure Technologies, Inc., all rights reserved.