
English – Instructions for sample collection2
Français – Instructions de prélèvement de l'échantillon3
Español – Instrucciones para la recogida de la muestra4

Contents: 1 kit containing 4 mL stabilizing liquid.

Summary and explanation of the kit:

OMNigene™-GUT DNA and RNA provides the materials and instructions for collecting and stabilizing microbial DNA and RNA from a fecal sample.

Warnings and precautions:

- FOR EXTERNAL USE ONLY.
- Do NOT remove the yellow tube top from the tube.
- Do NOT spill the stabilizing liquid in the tube.
- Wash with water if stabilizing liquid comes in contact with eyes or skin. Do NOT ingest. See SDS at www.dnagenotek.com.
- If collecting a liquid fecal sample, see separately provided user instructions.
- Small items may pose a choking hazard.

Storage: Room temperature storage pre- and post-collection 15°C - 30°C (59°F - 86°F).


Ship in accordance with applicable regulations covering transport of biological specimens.

Label legend:

- Collect sample by (Use by)
 - Catalog number
 - Lot number
 - Manufacturer
 - Caution, consult instructions for use
 - Consult package insert
- 15°C - 30°C Storage instructions
(59°F - 86°F)

ENGLISH

Proceed to transferring the fecal sample as soon as possible after defecation:


- 1 IMPORTANT PREPARATIONS:**
- Empty your bladder before beginning the collection.
 - Collect fecal sample free of urine or toilet water.
 - Utilize any toilet accessory that may be supplied with the collection device.
 - Toilet paper or tissues may be required.

- 2** While holding the yellow tube top, unscrew ONLY the purple cap and set aside for later use.
- IMPORTANT:** Do NOT remove the yellow tube top. Tube is not designed to stand upright. Do NOT spill the stabilizing liquid in the tube.


- 3** Use the spatula to collect a small amount of fecal sample.


- 4** Transfer the fecal sample into the top portion of the yellow tube top. Repeat until the sample fills the yellow tube top.
- IMPORTANT:** Do NOT push sample into the tube.

- 5** Scrape horizontally across the tube top to level the sample and remove any excess. Wipe exterior of tube and top with toilet paper or tissue as needed.


- 6** Pick up the purple cap and screw onto the yellow tube top until tightly closed.

- 7** For a **minimum of 30 seconds**, shake the sealed tube as hard and fast as possible in a back and forth motion.

- 8** Place spatula in original packaging or wrap in toilet paper and discard in garbage.
- IMPORTANT:** Send the sample for processing following the delivery instructions supplied separately by the kit provider.

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