

Optimizing high quality samples with OMNigene®•SPUTUM for *Mycobacterium tuberculosis* culture

Bitapi Ray, Jacques Niles, Andy Stewart and Cassandra Kelly-Cirino
DNA Genotek, Ottawa, Ontario, Canada

Challenge

Culture-based techniques remain the gold standard for diagnosing Mycobacterium tuberculosis (MTb) infection and for drug susceptibility testing. Sputum samples that are non-standardized, putrefied or otherwise low quality compromise culture results and impede patient diagnosis, and care. Conventional sputum decontamination procedures using NaOH/NALC can limit culture sensitivity and impact lab efficiency in several ways: (1) NaOH/NALC requires daily preparation and quality control, (2) the standard 15-minute processing can reduce viable MTb in samples by 20%-60%¹, (3) insufficient processing can lead to high levels of contamination, (4) excessive processing can significantly reduce MTb viability and lead to false negatives, and (5) processing-time constraints prevent batching of samples overnight or over weekends.

Solution

OMNigene®•SPUTUM provides a flexible alternative to NaOH/NALC processing and cold chain transport or storage. This simple-to-use reagent is a stable (1 year shelf-life), non-toxic chemistry that liquefies and decontaminates sputum while maintaining MTb viability (see results below). OMNigene•SPUTUM is added to samples at approximately 1:1 ratio and is used in two main ways: 1) at point-of-collection for reliable, easier transport of high-quality samples without cold chain, or 2) as a lab-added reagent that optimizes samples and allows for multi-day storage and batching without refrigeration.

In addition to stabilizing samples at ambient temperature, OMNigene•SPUTUM helps improve culture results by reducing putrefaction and culture contamination, and by curbing loss of viable MTb compared to that caused by NaOH/NALC treatment. Other important benefits of OMNigene•SPUTUM are that it simplifies laboratory workflows, integrates seamlessly with established diagnostic algorithms, and is compatible with solid and liquid (BBL® MGIT®) culture, as well as other MTb

diagnostic methods (e.g., smear microscopy, Cepheid® GeneXpert® MTB/RIF Assay, Hain Lifescience GenoType MTBDR^{plus}^{2,3}. Complementing these benefits, OMNigene•SPUTUM can help laboratories and national programs cut costs related to 1) cold chain stabilization during transport, and 2) the extra routine laboratory work that is required for NaOH/NALC preparation, processing and quality control.




Figure 1: Impact of OMNigene•SPUTUM on MTb viability

Method

Volumes of OMNigene•SPUTUM and phosphate-buffered saline (PBS), respectively, were spiked with either 105 cfu/mL (for 1-, 3- and 5-day time points) or 3×10⁶ cfu/mL (8-day time point) of attenuated MTb (aMTb; H37Ra) and were maintained at 4°C or 40°C for up to 8 days. A control solution (PBS+Tween) was also spiked as described and was stored at 4°C for 5 days. Each sample was centrifuged for 20 minutes at 3,500 × g to form sediment, which was re-suspended in 400 µL sterile PBS. From each re-suspended sediment preparation, 200 µL was inoculated into liquid culture, incubated at 37°C and checked daily. For the 1-, 3- and 5-day time points, cultures were grown in Middlebrook broth and viability was estimated visually against the McFarland turbidity standard. Cultures were called positive (i.e., “time to culture-positive” in days) when 0.5 McFarland was reached (approximately 1.5×10⁸ cfu/mL). For the 8-day time point, cultures were grown in BBL MGIT broth and were called positive when the BACTEC MicroMGIT reader triggered a positive reading.

Conclusion

OMNIgene•SPUTUM spiked with aMTb produces growth in liquid culture after 8 days at 4°C or 40°C. A minor delay in time to culture-positive (3-4 days) was observed in the OMNIgene•SPUTUM samples held at 40°C when compared to aMTB incubated in PBS.

Compatibility

OMNIgene•SPUTUM is compatible with liquid culture and maintains viable MTb for at least 8 days at 4°C to 40°C.

