



Primary packaging

## OG-500

### Features & Benefits

- Genomic DNA equivalent to DNA from blood for downstream applications
- Compatible with liquid handling and automated purification robots for high throughput
- Compact, robust design is ideal for storage and mailing
- Non-invasive, reliable self-collection that increases compliance and decreases costs
- Long term specimen stability at room temperature
- Higher quality and greater than 15x the DNA yield collected with buccal swabs
- Purification reagent available for easy purification
- FDA registered, CE certified

### Product Specifications

#### OG-500 (in packaging)

- Dimensions: 14.0 x 8.0 x 2.8cm
- Weight: 39g

#### OG-500 tube dimensions

- 16mm diameter
- 93mm height (without cap)



Instructions



Tube



Cap



Sample ready for transport/storage

