



New product brand to easily identify our DNA product from the RNA product that will be introduced soon.

Same PROVEN chemistry - ALL NEW Look.

We're always looking for ways to improve our product. Here's what's new based on customer feedback:



The easiest way to collect large amounts of DNA.




New plastics make the kit even easier to close and further protect the sample in transit to the lab.

New translucent base makes it easier for donors to see the "FILL TO" indicator. Lab technicians can now see the sample before opening the kit.

Upgraded user instructions: To maximize yield for every donor, we have removed the rinse step of the user instructions. This update eliminates any chance of a donor spitting the rinse water into the kit.

For more information or to download an editable version of the instructions, visit our website at:

http://www.dnagenotek.com/techsupport_donor.htm



Do NOT eat, drink, smoke or chew gum for 30 minutes before giving your saliva sample.

We'd love to hear from you. Send us your feedback at sales@dnagenotek.com