

Product change notification

Products affected: **ORagene-DNA** **ORagene-ONE**

Catalog number: OG-500, OG-510, OG-575 and ON-500

Description: Oragene®-DNA and Oragene-ONE for collection of human DNA

Effective date of change: October 2015

Customer action: Required if customer uses whitespace or barcode.

Change description:


Based upon user feedback and DNA Genotek’s continued commitment to quality, DNA Genotek has made improvements to our tube label. We will be aligning the date format with international standards, growing the label whitespace panel and repositioning label content to create additional writable whitespace. Additionally, the height of the machine readable barcode will be standardized to 7 mm and the font size on the human readable portion of the barcode will increase, making it easier to read.

	Existing	New
Left whitespace panel	15.3 mm × 38 mm	15.2 mm × 38 mm
Right whitespace panel	15.3 mm × 51 mm	20.8 mm × 60.1 mm
Date format	Month Year (e.g., Sept 2017)	YYYY-MM-DD (e.g., 2017-09-24)

Standard label

Notes:

Fill to line may vary depending on which Oragene kit is purchased. Images not to scale.






Standard label with barcode

Notes:

Fill to line may vary depending on which Oragene kit is purchased.

Images not to scale.



Reason for change:

Increase writeable whitespace based on customer feedback, change date format to align with international standards, and standardize machine and human readable barcodes for easy scanning and readability.

Please contact info@dnagenotek.com if you have any issues or concerns pertaining to this product change notification.