



# Product change notification

**Products affected:**  **ORAgene<sup>®</sup>DNA**       **ORAgene<sup>®</sup>DISCOVER**

**Catalog numbers:** OG-250 and OGR-250

**Description:** Oragene<sup>®</sup>DNA and Oragene<sup>®</sup>DISCOVER

**Effective date of change:** October 10, 2018

## Change description:

Oragene•DNA (OG-250) and Oragene•DISCOVER (OGR-250) disc format pre-collection storage instructions are changing from 15°C – 30°C to 15°C – 25°C. Additionally, the product shelf-life is changing from 30 months to 24 months from date of manufacture.

## Reason for change:

For optimal DNA, current data supports pre-collection storage temperatures between 15°C – 25°C for 24 months.

## Customer action:

Store pre-collected Oragene•DNA and Oragene•DISCOVER disc format kits between 15°C – 25°C for a 24 month pre-collection shelf-life.

Please contact [info@dnagenotek.com](mailto:info@dnagenotek.com) if you have any issues or concerns pertaining to this product change notification.