

# Product change notification

**Products affected:**  **ORAgene<sup>®</sup>Dx**

**Catalog numbers:** OGD-600, OGD-610, OGD-675  
OGD-500, OGD-510, and OGD-575

**Description:** Oragene<sup>®</sup>Dx

**Effective date of change:** This change will take place starting September 2019.  
Finished goods will be transitioned into the market after current inventory is depleted.

## **Change description:**

The product shelf life of Oragene-Dx saliva collection devices is increasing from 24 months to 30 months from date of manufacture.

## **Reason for change:**

Current data supports pre-collection storage temperatures between 15°C – 30°C for 30 months.

## **Customer action:**

Customers can now store pre-collected Oragene-Dx devices between 15°C – 30°C for a 30 month pre-collection shelf-life.

Please contact [info@dnagenotek.com](mailto:info@dnagenotek.com) if you have any issues or concerns pertaining to this product change notification.