



Product Change Notification

Product affected:

Part #: PG-100 and OA-400

Description: Performagene™•LIVESTOCK collection kit and Oragene®•ANIMAL collection kit

Effective month of change: May 2013

Customer action: Only required if PG-100 is integrated with automation systems.

Change description:

DNA Genotek is standardizing their animal products into one platform, therefore Oragene®•ANIMAL (OA-400) will be discontinued and the Performagene•LIVESTOCK (PG-100) name will change to Performagene. As well, the double-ended cap and the collection tube will change to the following dimensions:

Dimensions (mm)	New	Existing
Cap diameter	18.1 mm	20.1 mm
Collection tube height	92.1 mm	96.6 mm
Depth of tube, from top of tube to bottom of false-bottom	52.1 mm	52.3 mm
Depth of tube, from bottom of tube to bottom of false-bottom	40.0 mm	44.3 mm
Outside diameter (OD) of tube	15.0 mm	15.5 mm
Inside diameter (ID) of tube at threaded end	13.0 mm	13.4 mm
Inside diameter (ID) of tube at bottom of tube	13.5 mm	12.6 mm

Reason for change:

1. Based upon customer feedback, the utility of the PG-100 kit was easier to use with companion animals.
2. Improve lab workflow for those labs processing both PG-100 and OA-400 samples.
3. The PG-100 kits are less expensive than the OA-400 kits providing cost savings to our OA-400 customers.

Note:

We recommend that automated liquid handling protocols be updated to reflect the new tube dimensions, but expect that no other re-validation work will be required as a result of this enhancement.

Please contact info@dnagenotek.com if you have any issues or concerns pertaining to this Product Change Notification.