

Product change notification

Products affected: OMNIgene™•GUT, OMNImet™•GUT, PERFORMAbiome•GUT.

Catalog number/SKU: OMR-200, ME-200, PB-200, OM-200 and all customization with such SKUs

Effective date of change: Sept 2024

Change description:

To improve product supply to customers, better support workflow scalability, and improve sample traceability, DNA Genotek™ is changing suppliers for plastic cap product components. In addition, the new purple cap will no longer feature the DNA Genotek™ name. The change impacts OMNIgene™•GUT (OM(R)-200), OMNImet™•GUT (ME-200), AND PERFORMAbiome•GUT (PB-200). This change will be in effect for all products mentioned with a manufacture date after Sept 8th, 2024.

Current



New (in Sept 2024)



Product Weight Update

SKU	Pre-use device net weight/tare weight (g)	New pre-use device net weight/tare weight (g)	Pre-use with packaging (g)	New pre-use with packaging (g)
ME-200	13.77	11.99	28.8	22.41
OM-200	13.82/13.69	12.08	22.99	21.25
OMR-200	13.82/13.69	12.08	22.99	21.25
PB-200	13.82	12.08	22.0	20.26

Customer action: No action is required. Contact your DNA Genotek™ account manager or email info@dnagenotek.com if you have any issues or concerns pertaining to this product change notification.