

## Oral sample collection

**EASY • PAINLESS • PROVEN**

## Collect superior samples for your genetic analysis

The ORAcollection™RNA collection kit is an all-in-one system for the collection, stabilization and transportation of RNA from human saliva for Research Use Only. The user-friendly kit design makes it ideal for self-collection and assisted collection from donors of any age, either in a clinic or at home. The ORAcollection™RNA device stabilizes RNA at room temperature for months and is compatible with downstream applications such as RT-qPCR, microarray and miRNA sequencing.

### Applicable for:

- Human mRNA and miRNA expression profiling
- Developing "expression signatures" of disease
- Researching RNA found in saliva (i.e., human, viral and bacterial RNA)

### Specifications:

- Catalog number: OR-100
- Sample stability: 60 days
- Shelf life: 24 months



*Superior samples • Proven performance*

Patent ([www.dnagenotek.com/legalnotices](http://www.dnagenotek.com/legalnotices))

Some DNA Genotek™ products may not be available in all geographic regions.

ORAcollection is a trademark of DNA Genotek Inc., registered in various jurisdictions. All other brands and names contained herein are the property of their respective owners.

© 2024 DNA Genotek Inc., a subsidiary of OraSure Technologies, Inc., all rights reserved. PD-BR-364 v5/2024-09

**For Research Use Only**  
Not for use in diagnostic procedures

DNA Genotek Inc.  
3000 - 500 Palladium Drive  
Ottawa, ON, Canada K2V 1C2  
Subsidiary of OraSure Technologies, Inc.

Toll-free (North America): 1.866.813.6354  
Tel: +1.613.723.5757 • Fax: +1.613.723.5057  
[info@dnagenotek.com](mailto:info@dnagenotek.com)  
[www.dnagenotek.com](http://www.dnagenotek.com)

**DNAGENOTEK™**