

Product change notification



Product SKU: OCD-100, OCD-100A, OCR-100, OC-175, ORE-100, OR-100

Description: ORACollect•Dx, ORACollect•DNA, ORACollect for pediatrics, ORACollect•RNA

Effective date of change: Changes will be implemented in July 2020. Please note that it is possible to receive product after the July implementation date without the changes in effect. This may occur while existing inventory becomes depleted.

Change description:

The 1D 128C barcode found on the ORACollect tube label is a unique identifier which contains data related to the manufacturing of the device. The current barcode is comprised of 14 digits with the following barcode structure and meaning:

4	ORACollect product family
1X	Device product number
YYMMLL	Device lot number (year, month and lot number)
SSSSS	5 digit unique serial number

This product change notice is to inform of a change to the barcode structure. The new barcode structure will maintain a total of 14 digits, however the number of digits reserved for the device product number will change from 2 digits to 1 digit and the device lot number will change from 6 digits to 7 digits.

Current 14 digit barcode structure



New 14 digit barcode structure



Note: pink denotes variable content.

Reason for change: We are implementing this product change to increase manufacturing capacity to ensure product availability.

Customer action: No action is required.

Please contact info@dnagenotek.com or your Account Manager if you have any issues or concerns pertaining to this product change notification.