

# Product change notification

**Products affected:**  **OMNIgene<sup>®</sup>•ORAL**  **OMNIgene<sup>®</sup>•VAGINAL**  **OMNIgene<sup>®</sup>•SKIN**

**Catalog numbers/SKUs:** OMR-110, OMR-120, OMR-130 and OMR-140


**Description:** OMNIgene•ORAL (gum/plaque), OMNIgene•ORAL (tongue), OMNIgene•VAGINAL and OMNIgene•SKIN

**Effective date of change:** Early September 2021

**Change description:** The current white label on the OMNIgene collection tubes for each of the SKUs mentioned above will be replaced with a clear label and the barcode format will be updated to a unique 14 digit 2D barcode. In addition, the cap will be replaced with one without internal threads. Decapping/capping of the tube by the end user or laboratory personnel will not be hindered.

## OMR-110/120/130/140

Current label



Starting early September 2021



\*Pink text indicates variable information. The final label is black.

Starting early September 2021

14 digit 2D barcode -----



NNYYMMLLSSSS

2 digit product identifier code (11 for OMR-1XX SKUs) -----

YYMMLL: Year, month and lot number -----

SSSS: 5 digit unique serial number -----

**Reason for change:** DNA Genotek is updating manufacturing processes to improve product supply to customers, better support workflow scalability and improve sample traceability.

### DNA Genotek actions:

We are:

- moving manufacturing of these products to a new facility in order to implement the necessary changes and improvements to this product line;
- changing suppliers for certain product components to ensure better inventory management and integration into the new manufacturing processes.

As a result of these actions, production has been affected as we are experiencing some delays in implementing the new processes. We are working diligently to minimize any impact to our customers and apologize for any inconvenience.

**Customer action:** No action is required. Contact your DNA Genotek account manager or email [info@dnagenotek.com](mailto:info@dnagenotek.com) if you have any issues or concerns pertaining to this product change notification.