

# Product change notification

**Products affected:**  **ORAGENE-Dx**  **ORAGENE-DISCOVER**  **ORAGENE-DNA**  **ORAGENE-ONE**


**Catalog numbers/REFs:** OGD-600, OGD-610, OGD-675, OGR-600, OGR-675, OG-600, OG-610, OG-675, ON-600

**Effective date of change:** December 2022

## Change description:

A design change has been made to Oragene® 600 series cap. The cap wall thickness has been reduced from 3.7 mm wide to 1.3 mm. All other cap specifications remain unchanged.

## Oragene



## Reason for change:

The Oragene 600 series cap design has been implemented to support improved moulding capacity and product inventory.

## Customer action:

- No action is required for laboratories manually decapping collected Oragene 600 series samples.
- The cap design change is not expected to impact compatibility with automated decapping instruments; however, laboratories using automated decapping instruments should review the cap specifications to confirm there is no impact. To request caps to verify compatibility, contact your account manager.

Contact [info@dnagenotek.com](mailto:info@dnagenotek.com) if you have any issues or concerns pertaining to this product change notification.