

## Product change notification

**Products affected:**  OMNIGENE™•ORAL  OMNIGENE™•VAGINAL  OMNIGENE™•SKIN

**Catalog numbers/SKUs:** OMR-110, OMR-120, OMR-130 and OMR-140


**Effective date of change:** January 2024

### Change description:

As part of DNA Genotek's continuous improvement process, two changes are being made to the OMNIGENE™•ORAL, OMNIGENE™•SKIN and OMNIGENE™•VAGINAL product lines.


1. The OMNIGENE collection tubes have been available with a cap that does not have internal threads to grab the swab (known as the functional grabber feature). The updated OMNIGENE™•ORAL, OMNIGENE™•SKIN and OMNIGENE™•VAGINAL product line will now be manufactured with a cap that has the functional grabber feature for the swab. The collection tube has also been updated to work with the new grabber cap. Decapping/capping of the tube by users or laboratory personnel will not be hindered.
2. The current clear label on the OMNIGENE device collection tubes for each of the affected SKUs is being replaced with a white label for improved readability.

#### Current label



\*Pink text indicates variable information. the final label is black.

#### New label



\*Pink text indicates variable information. the final label is black.

### Reason for change:

DNA Genotek is updating manufacturing processes to implement improvements to this product line.

### Customer action:

There is no action required. All affected products manufactured after January 2024 will contain these changes. If you have any questions or want to discuss this product change notification, reach out to your account manager or email [info@dnagenotek.com](mailto:info@dnagenotek.com).